

Surface disinfectant for haemodialysis machines

Composition

100g of ClearSurf® contain:
9.2 g of Coconut-propylendiamine-guanidinium diacetate and 16 g of Benzalkoniumchloride

Properties

- Bactericidal and fungicidal;
- It has a rapid hepato-virucidal action. The inactivation of Hepatitis B and Hepatitis C viruses (surrogate BVD virus) is achieved after 30 minutes exposure time with a 0.5% solution or after 10 minutes with a 1% solution;
- It is odourless;
- Free from aldehydes or phenol and does not form aerosols;
- Suited also for cleansing and disinfection of PVC flooring, linoleum or acrylic.

Handling

ClearSurf is diluted in concentrations of 0.25%, 0.5% or 1%. The ready-to-use solution should be prepared on the same day it will be used.

Examples of concentrations and exposure times:

Bacteria and fungi: 1% - 15 min; 0.5% - 1 h; 0.25% - 4 h
Inactivation of HBV/HCV/HIV: 0.5% - 30 min; 1% - 10 min

Rota viruses Inactivation: 0.25% - 15 min

Presentation

6 x 2L bottles

