

DALI - ART

Therapy

LDL Aphaeresis

Flow rates

Blood flow: 10 – 500 ml/min

Citrate flow: 0 – 10 ml/min

Extracorporeal circuit protective systems

Monitoring of the arterial, venous and the adsorption column pressure

Air detector

Blood leak detector

Citrate drop detector

Removal rate

LDL Cholesterol: up to 76%

Lp(a) Lipoproteins: up to 70%

Cholesterol HDL: 8.6%

Blood/Citrate flow ratio

1:6 to 1:50

DALI columns volume / Total extracorporeal volume

500 mL (330 ml)

750 mL (430 ml)

1000 mL (580 ml)

Battery

Autonomy: approx. 15 minutes