

The Innovative material for PD

Our CAPD and APD systems are manufacture in Biofine®. Developed by Fresenius Medical Care, this innovative material is free from plasticizers and incinerates to harmless carbon dioxide and water.

It protects both the patient and the environment.

- Todos os componentes do sistema são produzidos em Biofine® o qual é isento de PVC e de DEHP.
- All system components are made from Biofine® which is free from DEHP and PVC.
- Advanced foil and bag design ensures durability of use.
- Lower water vapour loss during bag storage preserves the solutions better.
- Less adsorption and absorption of therapeutic agents, such as insulin, when added to the bag solution.
- Excellent clarity facilitates fluid bag inspection.
- Biofine® can be recycled.
- Biofine® is environmentally friendly and contains no chlorine, degrading simply to carbon dioxide and water when incinerated.

