

## Dialysis Fluid Purity

The “European Best Practice Guidelines” (EBPG) for haemodialysis recommend the use of pure water complying with the standards of the European Pharmacopoeia for the production of dialysis fluid. Starting with pure water, the ultra pure dialysis fluid is possible by using a dialysis fluid filter, DIASAFE®plus.

### DIASAFE® plus - Technical Data

#### Membrane material

Fresenius Polysulfone®

#### Surface area

2.2 m<sup>2</sup>

#### Housing Material

Polypropylene

#### Potting material

Polyurethane

#### Filter Dimensions

Length 28 cm/diameter 5.5 cm

#### Filtration rate

300 ml/h x mm Hg (3.75 l/min; max. 2 bar)

#### Endotoxin retention

≥ 106

#### Period of use / Filter exchange

max. 12 weeks or 100 ONLINE plus treatments or if integrity test fails

