

5008 Bloddlines



AV-Set ONLINEplus 5008
Product Number: 5018441

Arterial-Venous Set for Haemodialysis, ONLINE Haemodiafiltration and ONLINE Haemofiltration bipunction

Includes SafeLine for the substitution fluid priming necessary to carry out ONLINE Haemodiafiltration and ONLINE Haemofiltration treatments. It also allows for the ONLINE priming and cleaning of the extracorporeal circuit as well as the administration of ONLINE bolus and ONLINE blood reinfusion at the end of the treatment.

Filling volume: **132 ml**

Sterilisation method: Sterilisation by accelerated electrons (**e-beam**)

Units per box: **24 units**

AV-SN Set ONLINEplus 5008
Product Number: 5018511

Arterial-venous set for Haemodialysis, ONLINE Haemodiafiltration and ONLINE Haemofiltration **unipunction**.

Includes SafeLine for the substitution fluid priming necessary to carry out ONLINE Haemodiafiltration and ONLINE Haemofiltration treatments. It also allows for the ONLINE priming and cleaning of the extracorporeal circuit as well as the administration of ONLINE bolus and ONLINE blood reinfusion at the end of the treatment.

Filling volume: **166 ml**

Sterilisation method: Sterilisation by accelerated electrons (**e-beam**)

Units per box: **18 units**