

Prometheus[®]



Therapies

FPESA (Fractioned plasma separation and adsorption)

HD (Continuous dialysis)

HDi (Continuous dialysis)

SLEDD (Slow Dialysis)

Flow rates

Blood flow: 10 – 600 ml/min

Dialysis fluid flow: 0 – 300 – 500 – 800 ml/min

UF rate setting range: 0 – 4 L/h

Plasma flow: 10 – 500 ml/min

Extracorporeal blood circuit and protective systems

Monitoring of Arterial and Venous Pressure, PTM, Pre-Filter, Pos-Filter, Pse1, Pse2, Pads and TMPf

Air detector

Blood leak detector

Syringe Pump

Syringe delivery range: 0 – 10 ml/h

Flow: 0.1 ml/h

Maximum treatment duration

10 h

Battery

Autonomy approx. 15 minutes